

S3 CASTOR WHEEL ON PANTHERA X



Panthera S3 castor wheel




Panthera X with std castor


Information of what parts you need if you like to use the **new S3 castor wheel on Panthera X**.

These parts **MUST** be used when mounting the S3 castor wheel on Panthera X. Regardless if the wheelchair is on a new order or if you will change the castors on an existing wheelchair.


2100201
Gaffel 75, hjul 90, tapp X2
Fork 75, wheel 90, axle X2




2100204
Gaffel 90, hjul 90, tapp X2
Fork 90, wheel 90, axle X2



2100205
Gaffel 75, hjul 120, tapp X2
Fork 75, wheel 120, axle X2



2100206
Gaffel 90, hjul 120, tapp X2
Fork 90, wheel 120, axle X2



90 mm wheel

Chassis	Rear wheel	art.no.	Description
X	24"	2100201	Fork 75, wheel 90 mm, axle X2/24"
X	25"	2100201	Fork 75, wheel 90 mm, axle X2/24"
X	26"	2100204	Fork 90, wheel 90 mm, axle X2/24"

120mm wheel

Chassis	Rear wheel	art.no.	Description
X 5"	24"	2100205	Fork 75, wheel 120 mm, axle X2/24"
X	25"	2100205	Fork 75, wheel 120 mm, axle X2/24"
X	26"	2100205	Fork 75, wheel 120 mm, axle X2/24"
X 5"	26"	2100206	Fork 90, wheel 120 mm, axle X2/24"